

EDF-5005

Method 8290 Sample Fortification Solution

1.2 mL
in Nonane

Labeled Dioxins & Furans	Concentration (ng/mL)
2,3,7,8-TCDD (¹³ C ₁₂ ,99%)	100
2,3,7,8-TCDF (¹³ C ₁₂ ,99%)	100
1,2,3,7,8-PeCDD (¹³ C ₁₂ ,99%)	100
1,2,3,7,8-PeCDF (¹³ C ₁₂ ,99%)	100
1,2,3,4,7,8-HxCDF (¹³ C ₁₂ ,99%)	250
1,2,3,6,7,8-HxCDD (¹³ C ₁₂ ,99%)	250
1,2,3,4,6,7,8-HpCDD (¹³ C ₁₂ ,99%)	250
1,2,3,4,6,7,8-HpCDF (¹³ C ₁₂ ,99%)	250
OCDD (¹³ C ₁₂ ,99%)	500